Worldwide Acknowledged Skin Testing Equipment

Made in Germany



Multi Skin Test Center



Comprehensive Analysis of the Skin with Computer Support

Measurement of Skin Parameters

- World market leader for skin testing equipment since 1986
- Instruments used around the globe in research, dermatology and at the point of sale
- Full range of measurement parameters and cameras for skin type analysis
- From small handheld battery powered devices to smart multi-parametric systems with sophisticated software
- Ideal tools to promote cosmetic products & treatments
- Flexible customized solutions of hard- and software to suit the needs of your brand
- Combination of measurements according to your specific requirements
- From adding just your logo to creating your own design, everything is possible




Courage + Khazaka electronic GmbH

Mathias-Brüggen-Str. 91 50829 Köln Germany Phone +49 221 9564990 Fax +49 221 9564991 info@courage-khazaka.de www.courage-khazaka.de

Take Your Cosmetic Consultation to the Next Level with the Multi Skin Test Center and Skin Analysis Software SCS





Software Skin Competence System

- Wow effect of the vivid, meaningful & detailed analysis for your customers
- Facial Skin Tour or free **face & body** measurements
- Address the top skin concerns of your customers
- Individual recommendation of products
- Completely intuitive & full of possibilities
- All data, results & images can be saved
- Customize the program to your Cl (name, color, logo, images or videos)
- Also supporting touch screen devices



measurement moisture (cheek)



measurement sebum (cheek)




measurement elasticity (temple)



measurement redness (forehead)



measurement pigmentation (neck)



measurement barrier function (arm)

Multi Skin Test Center MC 1000

• multifunctional device

• moisture, sebum, pigmentation & redness, elasticity measurements

• Visioscope® PC 35 camera





Multi Skin Test Center MC 750


- multifunctional device
- model B2 with moisture, sebum,
 elasticity, pigmentation & redness measurement, further models with different probe configurations
- optional probes: temperature, pH, barrier function

Competent Skin Advice at the Point of Sale

Visioscope® PC 35

- USB-camera with parallel & cross polarized light to monitor skin surface and hair
- calculation of wrinkles, spots, pores, evenness of the complexion
- free image taking of the body
- with the special foils Sebufix® F 16 & Corneofix®F 20 for the determination of oil production and scaliness





Sebufix® & Corneofix®



determination of sebum production



evaluation of crow's feet

Visiopor® PP 34N



ra with special
UV light and polarisation filter

PORPHYRIN FLUGRESCENCE MEDIUM

PORPHYRIN FLUGRESCENCE MEDIUM

Minication for balanced kind confant

Minicati

determination of porphyrines

to monitor
fluorescent
porphyrines
(acne bacteria)



evaluation of spots/discolorations



evaluation of the visibility of pores