

## What does it measure?

The Meibometer® MB 560 is a unique tool to measure the **sebum content of the lacrimal fluid**.

## The Measuring Principle

The measurement is based on **grease spot photometry**. The Meibometer® strip consisting of Sebumeter® tape is brought into contact with the lacrimal fluid on the lower eyelid margin. It becomes transparent in relation to the lipid content of the tear film. Then the strip is inserted with the slider into the device and the **transparency of the strip**, where it had been in contact with the lacrimal fluid, is **measured by a photocell**. The peak of the light transmission curve represents the sebum content.

## Fields of Application

There are several applications in human and veterinary ophthalmology where the sebum content of the lacrimal fluid is of interest.

- It is important for **basic research** on the meibomian glands and the tear film.
- It is used for basic research of eye diseases, especially regarding the **dry eye syndrome**.
- It is important for **efficacy testing** of pharmaceutical products and **safety testing** of cosmetics used around the eye.


## Advantages

- Quick and **easy handling** of the foil strip for collecting the lacrimal fluid.
- Comfortable and **reproducible measurement** with the strip slider of the device.

- Display of the lipid content results in the **MPA CTplus software**.
- Interface and power supply by **USB connection**, no extra power supply needed.

## Additional Feature for the Scalp

- With an optional applicator the strips can be brought onto the **scalp with parted hair** to assess the **sebum content**.
- Ideal tool for **R&D of hair care products** (claim support and efficacy testing) as well as **dermatological basic research**.



## Technical Data

Dimensions: 13 x 5 x 18.2 cm (+ 11 cm with extended slider), Power supply: via USB, type B connector; Weight: 0.9 kg

Measurement principle: photometrical, Strip: Sebumeter® foil

Technical changes may be made without prior notice.

Courage+Khazaka electronic GmbH since 1986  
Mathias-Brüggen-Str. 91 · 50829 Köln · GERMANY  
phone +49 221 95 64 99 0 · fax +49 221 95 64 99 1  
info@courage-khazaka.de · www.courage-khazaka.de

